




Electromagnetic Healing.
Go Beyond the Surface.

Fast Forward Your Implants Practice


Magdent's innovative pulsed electromagnetic field (PEMF) healing abutments go beyond the surface to stimulate and reinforce the body's natural healing process for dental implant procedures.

The patented technology accelerates regeneration of bone and soft tissue, shortening osseointegration time, and enabling successful implant procedures in high-risk patients like heavy smokers, diabetes and osteoporosis patients, etc..

By reducing inflammation, the technology becomes an essential component for suppressing peri-implantitis and accelerating bone growth.



Dental Implantation Challenges




The Long Healing Period

This process can take between 3-6 months (and up to 9 months in cases of bone grafts)



Failure with High - Risk Patients

patients with osteoporosis, diabetes, and heavy smokers (~20%) have a high failure risk



Peri-implantitis

inflammation appears in around the dental implant in ~43% of cases, and can cause implant failure

There Is a Solution:

PEMF Bone Growth Stimulators (BGS) accelerates bone healing by up to x3

Pulsed Electro Magnetic Field (PEMF)—low intensity electric and magnetic fields combined. It is generated by passing an alternating current through a coil.

PEMF creates potential difference at the cell membrane level and stimulates the differentiation and proliferation of stem cells (MSCs) into bone forming cells (osteoblasts).

Since **1979**, PEMF devices have been **FDA approved** for bone-growth stimulation (Osteogenesis)

Companies like **Zimmer -Biomet®** and **Orthofix®** sell PEMF devices under **HCPCS code E0747**

Orthofix® PEMF devices annual sales are **\$200M**

MED[®] Bone Growth Stimulator



Similar design

Shaped like a healing abutment

No learning curve

Used according to standard healing abutment protocols

Made of

Titanium body incorporating a battery, coil and microelectronics

Variety

Can be compatibility with up to 85 % of implants brands

Clinical Value

A unique solution for the dentist

Duration

30 days of continuous operation upon activation

Envision the future of implantology




Our Product


The MED[©] enhances osteogenesis after implantation

How it works:

Our Product:
**The MED[®] stimulates osteogenesis
with inflammatory condition**



Acceleration of the
healing process



Increase in bone
implant contact



Increase in trabecular
bone volume density



The only treatment
for Perimplantitis

Fully Commercialized Product which can fits to up to 85% of the brands sold globally

The future of implantology is here

Strong Scientific Backing

Published in Peer Reviewed Journals:

- 2 clinical trials
- 1 animal trial
- 2 in-vitro trials

In Progress:

- 2 clinical trials – Osteogenesis
- 2 clinical trials – Peri - implantitis
- 1 clinical trial - Keratinized Gingiva






Randomized controlled clinical trial on 40 implants placed in 20 patients

Implant stability change from baseline in ISQ (Implant Stability Quotient)



TNF α



Effect of the Pulsed Electromagnetic Field (PEMF) on Dental Implants Stability: A Randomized Controlled Clinical Trial

Materials. 2020 Apr 3;13(7):1667

- ✔ Significant 13% increase in implant stability at test group Vs. 2% decrease in control over all time points
- ✔ 50% difference in TNF α raise in 4 weeks after implantation
- ✔ 25% difference in radiograph evaluation of marginal bone loss




Retrospective controlled clinical trial on 24 implants placed in 12 patients


Upper Jaw

Lower Jaw

Evaluation of the Maxillary implant stability quotient: MED vs. Passive healing cups



Evaluation of the Mandibular implant stability quotient: MED vs. Passive healing cups



Miniaturized Electromagnetic Device Abutment Improves Stability of Dental Implants

Journal of Craniofacial Surgery. 2019 Jun;30(4):1055-1057.


- ✔ Significant difference in stability rates between groups
- ✔ Similar effect at all bone types

The implant stability quotients were calculated at 15, 30 and 50 days of healing before conventional loading. Results are expressed as the percentage change from the baseline to the loading time point \pm SD.


The implant stability quotients were calculated at 15 and 30 days of healing before conventional loading. Results are expressed as the percentage change from the baseline to the loading time point \pm SD.

Controlled preclinical study on New Zealand rabbit tibia with microCT and histology

Micro CT



Histology




A new device for improving dental implant anchorage: a histological and micro-computed tomography study in rabbit

Clinical Oral Implants Research 2016 Aug;27(8):935-42.


- ✓ 48% increase in Bone to Implant Contact (%OI)
- ✓ 62% increase in trabecular bone volume density (BV/TV)
- ✓ 73% increase in connectivity density (Conn. D)
- ✓ 44% increase in the number of trabeculae (b.N)
- ✓ 32% decrease in trabecular spacing (b.Sp)

Controlled in-vitro study to assess osteoblasts formation in normal & inflammatory conditions

New Mature Bone Cell Formation



1 days 3 days 15 days 20 days



Pulsed electromagnetic fields increase osteogenesis commitment of MSCs via the mTOR pathway in TNF- α mediated inflammatory conditions: an-in vitro study


Scientific Reports (Nature publishing) 2018 Mar 23;8(1):5108 .

- ✔ Increase in MSC differentiation to osteoblasts
- ✔ Increase in MSC differentiation to osteoblasts in inflammatory conditions
- ✔ Enhanced mTOR signaling by increasing AKT, MAPP kinase, and Rraga values
- ✔ Increase in Runx^{*} osteopontin osteonectin osteocalcin collagen type I, wnt[†] foxO, ALP, BMP2 and BMP7 levels
- ✔ Lower levels of *IL1B and higher levels of IL10



Controlled in-vitro study to assess the influence of (PEMF) on bacterial biofilm colonization around implants

Changes in bacterial biofilm around implant



Mean bacterial counts (x10⁵) of the biofilms formed on the titanium surfaces of the Control group, without the electromagnetic healing cap being activated, and the Test group, with the electromagnetic healing cap activated on the first day and kept for the seven days of biofilm formation. The data were analyzed using the Mann-Whitney test (p < 0.05; p < 0.01, p < 0.001).

Antimicrobial effects of a pulsed electromagnetic field: an in vitro polymicrobial periodontal subgingival biofilm model


Biofouling. 2020:1 -8

- ✔ The mean total bacterial counts were lower in the Test group vs the control group (p < 0.05)
- ✔ Antimicrobial effects on the bacterial species and can be used to control bacterial colonization around dental implants

Peri-implantitis Clinical Trial

A prospective double blind randomized controlled trial that examined PEMF in treatment of periimplantitis on patients that were treated with dental implants and crowns a few years ago

One month after treating with MED, we saw a reduction of marginal bone loss which remains constant after 3 months follow up



The clinical trial was completed during September 2022, and is being finalized for publication these days.

Focused Pulsed Electromagnetic Field is an Effective Adjunctive Treatment in implant with periimplantitis: prospective randomized double-blind controlled clinical trial

Yaniv Mayer, Juan Khouri, Jacob Horwitz, Hadar Zigdon

Department of Periodontology, School of Graduate Dentistry, Rambam Health Care Campus (RHCC), Haifa, Israel.

Introduction: As dental implants becoming the first choice to replace teeth that were lost or congenitally missing, an ever-greater concern for periimplantitis and its consequences is imminent. Potential use of PEMFs as modulator of immune responses alone or in combination with pharmacological therapies represents a novel frontier of investigation with interesting clinical perspectives.

The Miniaturized Electromagnetic Device (MED – which was invented, designed and owned by Magdent®), is small enough to fit most dental implants. The MED is shaped like a simple healing abutment and is easily screwed into the implant. It is removed at the end of the treatment period following standard healing abutment placement protocols. *The purpose* of this randomized, double-blind clinical trial was to compare traditional non-surgical implants' surface debridement alone, to non-surgical implants' surface debridement combined with pulsed electromagnetic field (PEMF) generation around dental implants with periimplantitis.

Electromagnetic Healing: Go Deeper, Heal Faster



Higher success rate with high risk patients
(heavy smokers, diabetes, oncologic patients, sponges bone etc.) - increasing patients pool



improves the **quality of care**, increase efficiency & **productivity** to your practice.



Offering **cutting-edge technology** without investing in capital equipment



Increase revenue by charging the patient extra for a **premium treatment method**



Stay on-top of the trend and technology
Improve your clinic reputation ahead of the curve


Accelerate Your Implants Procedure

Faster recovery


Successful implantation in high-risk patients

Reducing inflammation

Bone Growth in Inflammatory Condition



Age: 60 Years old; Patient Condition: Healthy; Procedure: Implants & Stable Bridge




18/7/2019
Insertion Day

16/09/2019
Impression Day


7/10/2019
1 Week post Crown restoration

61 days

20 days




Age: 71 Years old; Patient Condition: Healthy; Procedure: Implants & Stable Bridge




11/7/19
Extraction & Insertion Day

10/09/2019
Impression Day

60 days



Age: 63 Years old; Patient Condition: Healthy; Procedure: Implants & Stable Bridge




20/11/2019
Insertion Day

20/01/2020
Impression Day

30/01/2020
Restoration day

60 days

15 days



Smokers



Meta-Analysis > J Clin Periodontol. 2020 Apr;47(4):518-528. doi: 10.1111/jcpe.13257.


Epub 2020 Feb 7.

Levels of smoking and dental implants failure: A systematic review and meta-analysis

Roohollah Naseri¹, Jaber Yaghini², Awat Feizi³

Affiliations + expand

PMID: 31955453 DOI: 10.1111/jcpe.13257



Results: Having additional information supplied by the authors, 23 articles were selected for final analysis. The meta-analyses based on implant- and patient-related data showed a significant increase in the RR of implant failure in patients who smoked >20 cigarettes per day compared with non-smokers (implant based: $p = .001$; RR: 2.45; CI: 1.42-4.22 and patient based: $p < .001$; RR: 4; CI: 2.72-5.89).

Conclusion: The risk of implant failure was elevated with an increase in the number of cigarettes smoked per day.

Age: 70 Years old; Patient Condition: Smokers; Procedure: Implants & Stable Bridge




Heavy smoker, 6 implants in edentulous mandible. Magdent allowed a faster healing process


29/10/2019
Insertion Day

15/12/2019
Impression Day

45 days



Low bone density around implant



Age: 55 Years old; Patient Condition: Smokers; Procedure: Implants & Stable Bridge



Heavy smoker, close sinus lifting, including ridge preservation, a complete healing and final restoration within 82 days.

15/10/2019

3 Months after extraction

07/11/2019



Implantation & Crestal & Sinus Lifting

27/01/2020

Final restoration

24 days

82 days



Age: 73 Years old; Patient Condition: Diabetes;



Long term diabetes patient, Glucose 200 mg/dL, all teeth at lower jaw were extracted; D4 Bone, Immediate insertion of 7 implants



07/08/2019
Extraction & Insertion Day

24/10/2020
Impression Day

07/11/2019
Final restoration

75 days

14 days



Age: 50 Years old; Patient Condition: Patient suffers from Cancer, received biological treatment






Extraction & Insertion in the same day, 4 implants to support denture

26/09/2019
insertion day

11/11/2019
Impression Day

18/11/2019
Final restoration
7 days

46 days



Age: 55 & 70 Years old; Patient Condition: Peri Implantitis;



Patients suffered from peri implantitis, approximately 2 years post implantation. MED was used as an adjunctive treatment and witnessed bone growth even during inflammation around the implant.

Case 1 - 30 days





Case 2 - 90 days

Baseline

1 month

Baseline

3 month



Regulatory Approvals

CE


Annex II without section 4

This certification is subject to surveillance by MEDCERT.

Effective date: 2020-05-22
Expiry date: 2024-05-27

Report No.: 7162FS09F
Process No.: QS - 7162
Certificate No.: 7162GB410200522

Hamburg, 2020-05-22


MEDCERT Certification Body
(Dr. Andreas Schich)

ISO

EN ISO 13485:2016

This certification is subject to surveillance by MEDCERT.

Effective date: 2020-05-22
Expiry date: 2023-02-26

Report No.: 7162FS09F
Procedure No.: QS - 7162
Certificate No.: 7162GB445200522

Hamburg, 2020-05-22


MEDCERT Certification Body

FDA



In process

Clinical Recommendations



Enhanced osteogenesis: Producing Pulsed Electro Magnetic Field (PEMF), which has a major effect on osteoblasts, resulting in bone formation even in inflammatory conditions, Enhancing osteogenesis.



Magdent 's MED is the only real solution for Diabetic patients/ Smokers/ Cancer patients / Implants failure sites as it reduce their procedure risk.



Low Magnetic Current: The hollowed body of the healing cup contains electronic components. The Pulsed Electro Magnetic Field generates low level current in the bone around the implant surface of up to 2 mm. It is most effective around the top 7 mm of the implant & soft tissue.



Temporary Patial Denture: Maximizing MED's procedure success and eliminating temporary rehabilitation leads to higher success rates than prolonged temporary rehabilitation



Identical Standard Protocol

- MED should only be tightened gently by hand until it can no longer be screwed on. Don't damage the MED.
- Not suitable for immediate loading.
- Use the MED during the rehabilitation process to maintain the gums' health.
- Don't cover the MED under the gingiva.
- The implant must be positioned at the bone level (crestal).
- Make sure there's no occlusal interference with opposing teeth.

Meet the team:



Prof. Shlomo Barak
D.M.D

Founder &
Chairman of the Board



Adv. Benny Barak
LLB., M.B.A

Chief Executive Officer



Elinor Zohar

Director of R&D and
Manufacturing



Ran Tuttnauer

Board Member



Dr. Oleg Dolkart
PhD Medical Sc.

Chief Medical Officer



Elad Yakobson
B.Sc. Ind. Eng.

Chief Operations Officer



Galit Lass
Economics

Director of Marketing and
Business Development



Tal Shahar
B.Sc. Elec. Eng.

Head of Product

Advisory board




Dr. Ophir Fromovich D.M.D.

Founder and owner of Alpha Tec, which was
sold to Nobel BioCare for \$5M in 2008



Prof. Jamil Shibli PhD

Vice Dean of Dental Research, University of
Guarulhos. Reviewer of the journals: International
Journal of Oral and Maxillofacial Implants, COIR,
Journal of Clinical Periodontology, Journal of
Periodontology and Clinical Implant Dentistry.



Healing smiles through innovation!



Официальный дистрибьютор в России

8 (800) 777-09-72

+7 (918) 977-12-12

+7 (918) 483-39-03

www.paltop.ru

manager@paltop.ru

drpmp@mail.ru